



Mshot

ML11 Features:

- Exquisite design, field of view is wide and clear.
- 3W LED, Pre-set lamp center, brightness continues adjustable.
- Entry level microscope
- Solar power bank is optional for outdoor using.

Specification

Item	Specification
Eyepiece	WF10X/18 (WF10X/20 is optional)
Head	30°inclination binocular,P.D. 50mm - 75mm,single type diopter adjustable 30°inclination trinocular,P.D. 50mm-75mm,single diopter adjustable,fix light division R:T:80:20(optional)
Objectives	Achromatic objective 4X/0.10, W.D.: 37.5mm (Infinity plan achromatic objectives are optional) Achromatic objective 10X/0.25,W.D: 6.55mm Achromatic objective 40X/0.65,W.D: 0.669mm Achromatic objective 100X/1.25, W.D: 0.198mm
Nosepiece	Quadruple nosepiece with inward tilt
Focusing	155mmx142mm double layer mechanical, moving range 76mmX50mm , accuracy 0.1mm
Stage	Low position hand wheel coaxial coarse and fine adjustment with limit stopper and lock device, stage height movement (coarse movement stroke: 25mm),fine focusing accuracy 0.002mm
Transmitted illumination	3W LED, Pre-set lamp center, brightness continues adjustable Abbey condenser,N.A.1.25, with changeable aperture diaphragm Filters: blue, yellow and green
Power line	External set wide voltage power adapter, input 100V-240V-AC50/60HZ, outputDC12V2A

Guangzhou Micro-shot Optical Technology Co. Ltd

Add: Room 506 Vanke Cloud Building A / 1933 Huaguan Road Guangzhou / China
WEB:www.m-shot.com/ sales@m-shot.com / TEL:0086-20-38262481



Mshot

ML11

Biological Microscope

