



MZ101

MshOt

Stereo microscope MZ101 is outstanding with super long working distance and wide zoom range, satisfy different zoom times requirement and easy to operate, ergonomics design provide comfortable use experience. Both transmitted light and reflected light are white LED light source.

FEATURES:

- Wide zoom range, optional with adding assistant objectives to enlarge zoom times.
- High brightness white LED lighting.
- Stable and smooth focus, light and handy.

Specification

Item	Specification
Eyepiece	Wide field WF10X/22mm
	Wide angle
Head	45° inclination, diopter adjustable
Objective	Continuously zoom adjustable
	Zoom range 0.7~6.3
	Zoom ratio 1:9
	Working distance 110mm
Field of view	Φ33.8~Φ3.49
Zoom times	7X~63X
Lighting	Reflected : 3W high brightness LED
	Transmitted : 3W high brightness LED
Base and focus	White light LED , brightness adjustable
	Push-pull phase contrast unit, working distance 55mm
	Green filter

MshOt

MZ101

Stereo Microscope

