


LCD Display Tablet Camera



Features:

- Windows 10 OS.
- Hardware:64GB, memory:2GB.
- Scientific CMOS sensor.
- Sensor size 1/2.3', effective resolution 16.0MP.
- Digital field of view 90%, convenience to imaging.
- Power adapter stable working.
- Come with image processing software.
- Auto and manual exposure and white balance.



Windows tablet screen camera	Item	MD28-W
	Screen size	11.6 inch
	Sensor size	1/2.3 inch
	Resolution	16.0 megapixels
	Frame rate	30fps@1080P
	Output port	USB2.0, HDMI, TF (Inset WIFI)
	Pixel size	1.335um×1.335um
	Power	DC-12V/2A
	UI	Digitization, Mouse controllable, User friendly design
Image analysis function	Marking	Point coordinates, cross lines, coordinate systems, text annotations
	Measurement	Length: line, fold line and curve. Distance: parallel line, point to line
	Circle	Radius set circle, two-point set circle, three-point set circle
	Geometric	Line length, radius, two/three points circle, concentric circle
	Geometric area	Polygon, Square
	Image snapshot	Support 16 megapixels image capture
	Record video	Support 1080P video