


Mshot

Mshot

Microscope Digital Camera

MD80-HD





MD80-HD

### Features:


- Ultra HD HDMI Camera
- 1400TV lines of resolution
- Wireless IR remote control

### Parameter :

|               |                        |
|---------------|------------------------|
| Item          | MD80-HD                |
| Sensor Type   | Sony Sensor 1/1.9 inch |
| Resolution    | 3840 × 2160            |
| Pixels Size   | 1.85×1.85um            |
| Frame Rate    | 30 FPS                 |
| Definition    | 4K UHD                 |
| Pixel depth   | 12bit                  |
| White Balance | Manual/ AUTO           |
| Exposure      | Manual/ AUTO           |
| Freeze        | Support                |
| Gain Control  | Support                |
| Photo & Video | Support                |
| Wide Dynamic  | Support                |
| Accessories   | IR remote control      |
| Power Supply  | 12V 2A                 |



Intelligent Imaging for Great Discoveries



4K UHD Microscope Camera



Guangzhou Micro-shot Optical Technology Co. Ltd

ADD: Rm 506, Lobby A, Vanke Cloud, 1933 Huaguan Road, Tianhe District, Guangzhou, 510640 China  
 WEB : <http://www.m-shot.com> Email: [sales@m-shot.com](mailto:sales@m-shot.com) TEL: 020-38250606/38262481

