

Touch screen MG-100

- ▶ High sensitivity touch screen design
- ▶ Know current band and brightness and adjusted according to needs
- ▶ Show light box temperature and working time
- ▶ Human body sensing function is optional, light off when people walk away




MG-30




- ▶ Modular and elegant design
- ▶ With high luminous efficiency LED white light beads, the luminous flux of 30W meets majority customer needs
- ▶ High linear light control device is flexible to experiments
- ▶ Strong applicability used as bright field / phase contrast / dark field and also BG fluorescence excitation

Specification

Model	MG-100	MG-30
Wavelength	350-390nm/410-760nm	410-760nm
Host power line	1.5 meter, 8P aviation head (With temperature control)	
Power line	According to customer needs	
Connector	Olympus, Nikon, Leica, Zeiss	
灯珠光谱图		

MshOt 明美

Broad-spectrum LED Light Source MG-100/MG-30 Digital Touch Screen



Guangzhou Micro-shot Technology Co.,Ltd

Add: Rm A-506, Vanke Cloud, 1933 Huaguan Road, Tianhe District, Guangzhou, China
Website: www.m-shot.com Email: sales@mshot.com Tel: 0086-20-38262481

Microscope Imaging Solution Provider